HEALTHCARE™

Clorox Healthcare™ Bleach Germicidal Cleaner

One-step cleaning, disinfecting and deodorizing of hard, nonporous surfaces*

For more information,
contact your Clorox sales
representative or call
800-234-7700.

email:
healthcare@clorox.com
visit us:
www.cloroxhealthcare.com

1. REMOVE gross soil if visible. For *C. difficile* spores, always clean surface prior to disinfecting.†

2. SPRAY 6"-8" from surface until surface is completely wet.

3. To DISINFECT, let stand for 1 MINUTE. To kill TB, allow 30 seconds contact at room temperature. To kill *C. difficile* spores, allow 5 minutes contact. WIPE with a clean, damp cloth. Allow to air dry.

This product may be used on stainless steel, glazed ceramic tile, hard plastics, nonporous vinyl, painted surfaces, Plexiglas®, laminated plastic countertops, fiberglass (sealed) surfaces, glazed porcelain and glass.

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling. Do not use this product with ammonia or acids such as vinegar, rust removers or toilet bowl cleaners. Not intended for use on fabrics. Do not use this product with ammonia or acids such as vinegar, rust removers or toilet bowl cleaners. Not intended for use on fabrics. Use by expiration date.

© 2012 Clorox Professional
Products Company,
1221 Broadway,
Oakland, CA 94612.

NI-18559

* Look on the product label for additional directions for other applications such as for food-contact surfaces, veterinary sites and instrument precleaning or disinfection. For *C. difficile* spores, EPA requires surface cleaning prior to disinfecting.

† This is a one-step product. Precleaning is required only if blood, vomit, fecal matter or other soil is visible, or in advance of disinfecting *C. difficile* spores.

Use product as directed.

Plexiglas is a trademark of Arkema France Corporation.

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© 2012 Clorox Professional Products Company,
1221 Broadway,
Oakland, CA 94612.

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Plexiglas es una marca registrada de Arkema France Corporation.